

Supplement II to JP17 table of errata

April 28, 2020

General Test / 2.51 Conductivity Measurement

Page	Line	Correction	Error
2891	left ↑ 1	$\kappa = 1000 \sum_i^{all\ ions} C_i \cdot \lambda_i$	$\kappa = 100 \sum_i^{all\ ions} C_i \cdot \lambda_i$

General Test / 9.21 Standard Solutions for Volumetric Analysis

Cerium (IV) Sulfate, 0.1 mol/L

Page	Line	Correction	Error
2901	left ↓ 28-30	<p><i>Preparation</i>—Dissolve 40.43 g of cerium sulfate (IV) tetrahydrate <u>in 550 mL of diluted sulfuric acid (1 in 11). After cooling, add water</u> to make 1000 mL, and standardize the solution as follows:</p>	<p><i>Preparation</i>—Dissolve 40.43 g of cerium sulfate (IV) tetrahydrate <u>in water</u> to make 1000 mL, and standardize the solution as follows:</p>