




Pyrantel Pamoate

A solution prepared as follows: Dissolve 0.1 g of Pyrantel Pamoate in 50 mL of dimethylformamide, and add methanol to make 200 mL. To 2 mL of this solution add the solution of hydrochloric acid in methanol (9 in 1000) to make 100 mL.

Pyrazinamide

A solution in 0.1 mol/L hydrochloric acid TS (1 in 100,000)

Pyridostigmine Bromide

A solution in 0.1 mol/L hydrochloric acid TS (1 in 30,000)